**de – Gebrauchsanweisung Ultraschallvernebler
en - Instructions for use Ultrasonic nebulizer**

USV 2010

1 Inhalt

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3.1	General safety information	19
3.2	Strictly follow these instructions for use	19
3.3	Maintenance	19
3.4	Accessories	19
3.5	Product-specific safety information	19
3.6	Patient safety	19
3.7	Information on electromagnetic compatibility.....	19
3.8	Inspection of delivery	20
3.9	Device symbols.....	20
4	Safety instructions	20
5	Device description.....	20
5.1	Device description	20
5.2	Intended use	20
5.3	Area of application and operating conditions	21
5.4	Self-evident improper use	21
5.5	Device components.....	21
5.6	Functional description.....	22
6	Putting into operation.....	23
6.1	Preparation	23
6.2	Connecting to power supply.....	24
6.3	Functional test.....	24
6.4	Preparing USV 2010.....	24
6.4.1	Nebulizing directly from nebulizer chamber	24
6.4.2	Nebulizing with standard disposable sterilized water systems.....	25
6.4.3	Nebulization from the disposable medication container.....	25
6.4.4	Nebulization with refilling system (optional).....	26
6.5	Assembling accessories	26
6.5.1	Mounting the table model (Optional accessories).....	26
6.5.2	Mounting the wall-mounted model accessory set (Optional accessories)	27
6.5.3	Mounting the trolley accessory set (Optional accessories)	27
6.5.4	Transport with trolley.....	28
7	Operation	28
7.1	Functional test.....	28
7.2	Starting operation	28
7.3	Stopping operation.....	28
8	Faults, troubleshooting	29
8.1	General.....	29
8.2	Replacing the quartz oscillator	29
8.3	Replacing the fuses	29
9	Reprocessing	30
9.1	Disassembly.....	30
9.2	Notes on implementation	30
9.3	Implementation.....	30
9.4	Manual cleaning of the quartz	31
9.5	Cleaning the automatic feeder set	31
10	Preventive maintenance / repair	31
10.1	General.....	31
10.2	Preventive maintenance	31
10.3	Repair	31
11	Available products	31
11.1	Standard delivery	31
11.2	Optional accessories.....	31
11.3	Spare parts and consumables.....	32
11.4	Interface description	32
12	Disposal.....	32
13	Technical data	32
13.1	Ambient conditions	32
13.2	Performance data.....	32
14	Electromagnetic compatibility	33
14.1	Electromagnetic emission	33
14.2	Electromagnetic interferences in other devices	33
14.3	Electromagnetic immunity	33
14.4	Electromagnetic immunity, devices not intended for life support	34
14.5	Recommended separation distances	34

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1 Safety information definitions

Warning

A WARNING statement provides important information about a potentially hazardous situation which, if not avoided, could result in death or serious injury.

Caution

A CAUTION statement provides important information about a potentially hazardous situation which, if not avoided, may result in minor or moderate injury to the user or patient or in damage to the medical device or other property.

Note

Provides additional information.

2 Definition of target groups

For this product, users and service personnel are defined as target groups. These target groups must have the necessary training and knowledge to use, install, reprocess, or maintain the product. The product must be used, installed, reprocessed, or maintained exclusively by defined target groups.

2.1 Users

Users are persons who use the product in accordance with its intended use.

2.2 Service personnel

Service personnel are persons who are responsible for the maintenance of the product. Service personnel must be trained in the maintenance of medical devices and install, reprocess, and maintain the product.

3 For your safety and that of your patients

3.1 General safety information

The following **WARNING** or **CAUTION** statements apply to general operation of the medical device. **WARNING** or **CAUTION** statements specific to subsystems or particular features of the medical device appear in the respective sections of these instructions for use or in the instructions for use of another product being used with this device.

3.2 Strictly follow these instructions for use

Warning

Any use of the medical device requires full understanding and strict observation of all sections of these instructions for use. The medical device must only be used for the purpose specified under Intended use on page 20. Strictly observe all **WARNING** or **CAUTION** statements throughout these instructions for use and all statements on medical device labels. Failure to observe these safety information statements constitutes a use of the medical device that is inconsistent with its intended use.

3.3 Maintenance

Caution

Risk of medical device failure and of patient injury.

The medical device must be inspected and serviced regularly by service personnel. If the above is not complied with, medical device failure and patient injury may occur. Observe chapter Preventive maintenance / repair.

3.4 Accessories

The manufacturer only assumes responsibility for the safety, reliability and usability of the device if original devices or original spare parts are used and only authorized persons repair the device. The statutory warranty period shall apply. In case of a warranty claim, please contact your authorized specialist dealer.

Warning

Only the accessories indicated on the order list have been tested and approved for use with the medical device. Therefore, it is strongly recommended that only these accessories are used in conjunction with the medical device. Otherwise, correct functionality of the medical device can be jeopardized and no warranty is given.

3.5 Product-specific safety information

Warning, risk of malfunction

Impermissible modifications to the medical device lead to malfunction.

This medical device may not be modified without permission from its manufacturer.

3.6 Patient safety

The design of the medical device, the accompanying documentation, and the labeling on the medical device are based on the assumption that the purchase and the use of the medical device are restricted to professionals, and that certain inherent characteristics of the medical device are known to the user. Instructions and **WARNING** and **CAUTION** statements are therefore largely limited to the specifics of the medical device. These instructions for use do not contain references to various hazards which are obvious to professionals who operate this medical device as well as references to the consequences of medical device misuse, and to potentially adverse effects in patients with different underlying diseases. Medical device modification or misuse can be dangerous.

3.7 Information on electromagnetic compatibility

General information on electromagnetic compatibility (EMC) according to international EMC standard IEC 60601-12.

Medical electrical equipment is subject to special precautionary measures concerning electromagnetic compatibility (EMC) and must be installed and put into operation in accordance with the EMC information provided on page 33.





Portable and mobile radio frequency communication equipment can affect medical electrical equipment.

3.8 Inspection of delivery

Upon delivery of the device, check if all components are included in the delivery and in good condition. Take the device out of its packaging and start by checking the delivery for completeness according to the order (11.1 Standard delivery). In case of customer claims, please contact your authorized specialist dealer.

3.9 Device symbols

The following symbols are used on the device:

	Follow the instructions for use
	Material name of the plastic
	CE symbol and number of the notified body according to the EC Directive 93/42/EEC for medical devices
	Off On
	Device with applied part, type B
	Alternating voltage
	Drip-proof
	Device fuse
	Marked in accordance with EC Directive 2002/96/EC (Waste Electrical and Electronic Equipment Directive) Symbol for "This device may not be disposed of at municipal collection points for waste electrical and electronic equipment."

4 Safety instructions

Warning!

The ventilation slots must not be blocked by objects, otherwise the device does not work properly and the aerosol nebulization will be disturbed.

Warning, risk of infection!

In order to avoid an infection or bacterial contamination, the relevant hygiene rules must be complied with. Observe the intended use of the bacterial filter. Always wear gloves when performing work.

Warning, risk of injury!

The heatable hose connection at the rear of the device and the patient must not be touched at the same time.

Damages

The device must be checked for damage (housing, etc.) and safety-related faults before each use. Damaged devices must be repaired by an authorized body before putting them back into operation. All works that require the use of tools must be performed by the manufacturer or bodies authorized by it.

Power supply

The mains voltage must comply with the information on the rating plate. The mains plug must be disconnected before cleaning. No fluids must penetrate into the inside of the device.

Liability

The manufacturer is not liable if the device is used for purposes other than those described in the instructions for use or if non-original spare parts are used without authorization, nor can any warranty claims be made. Only the manufacturer or authorized professionals are authorized to perform repair work, using original spare parts and complying with the applicable safety measures for medical electrical equipment.

Warning, observe hygiene regulations!

Contaminated device parts can endanger the health of patients. Before commissioning, the device must be reprocessed in accordance with the operational hygiene regulations.

Caution, functional inspection!

The user must check the functionality and suitability of the applied parts for the respective intended use as well as their biocompatibility.

5 Device description

5.1 Device description

The device consists of a washable plastic housing with interior electrical components and an attachable nebulizer chamber. It is commissioned and controlled with push buttons.

The foot of the nebulizer chamber contains a quartz oscillator that causes the fluid in the container to nebulize due to its high-frequency oscillations. A fan located in the housing moves the aerosol through the patient hose and into the desired respiratory area of the patient. Prior to nebulization, the aspirated air is led through a bacterial filter.

There are three different kits for the different ultrasonic nebulizer installation options:

- Table version: The device is located on a table or a level and robust surface
- Wall-mounted version: The device is mounted on a standard rail 25x10 mm (DIN)
- Trolley-mounted version: The device is mounted to a trolley pole

5.2 Intended use

The ultrasonic nebulizer USV 2010 is only intended for humidifying the breathing air of spontaneously breathing patients and for aerosol therapy. An aerosol therapy must only be carried out as directed by a medical professional.

The ultrasonic nebulizer USV 2010 is used at hospitals, medical practices, nursing homes, care homes and in the military medical service by specifically trained specialist personnel.

Aerosol therapy must only be carried out as directed by a medical professional

5.3 Area of application and operating conditions

The medical device may be used in medical areas.

Warning, risk of explosion

This medical device is neither approved nor certified for use in areas where the oxygen concentration may exceed 25 Vol% or combustible or explosive gas mixtures are likely to occur.

Warning, risk of explosion

Do not use the device together with combustible or explosive liquids.

Warning, risk of patient injury

Do not use the ultrasonic nebulizer on patients requiring ventilation

Warning, risk of patient injury

Do not install the ultrasonic nebulizer in breathing systems

Warning, risk of patient injury

Do not use the ultrasonic nebulizer with a ventilation mask, tracheal cannula, or mouthpiece. This can lead to accumulations of liquid in the respiratory system.

Warning, risk of infection

Do not use the ultrasonic nebulizer without a bacterial filter

The air flow is only filtered adequately if a suitable bacterial filter is used.

Replace the bacterial filter after 24 hours of operation or in accordance with the special instructions provided by the respective filter manufacturer.

Reusable hoses must be cleaned once a day or whenever there is a change of patient.

Disposable hoses must be replaced once a day or whenever there is a change of patient.

The ultrasonic nebulizer must only be used for its intended purpose. Due to the respirable particle size of the mist, particular care regarding hygiene must be taken during preparation and execution of the inhalation. Medication should only be nebulized with a medication container (11.2). A medical professional must determine the need for nebulizing, especially regarding the inhalation of medication.

5.4 Self-evident improper use


Warning, observe environment of use!

The device must not:

- be covered, thrown, forcefully opened or dropped,
- be used for patient ventilation,
- be used without bacterial filter or fluids,
- be operated by children who are unsupervised,
- be used together with combustible or explosive liquids or operated in potentially explosive areas.

No objects must be placed in the housing

5.5 Device components

Similar image

- | | | | |
|---|-----------------------------------|----|-----------------------------|
| 1 | Ultrasonic nebulizer | 7 | Bacterial filter connection |
| 2 | Nebulizer chamber, complete | 17 | Nebulizer chamber |
| 3 | Hose approx. 300 mm | 18 | Nebulizer chamber lid |
| 4 | Hose approx. 1,200 mm (accessory) | 19 | Automatic feeder set |
| 5 | Bacterial filter | 20 | Quartz module |


- | | | | |
|----|---|----|------------------------------|
| 12 | Connector for heatable hose (accessory) | 8 | ON/OFF push button |
| 13 | ON/OFF switch | 21 | Mist volume push button |
| 14 | Device socket | 9 | Air volume push button |
| 15 | Voltage selector | 10 | Timer push button |
| | | 11 | Holder for nebulizer chamber |
| | | 22 | Support |

5.6 Functional description

In the ultrasonic nebulizer USV 2010, a quartz oscillator generates high frequencies, which transform the fluids in the nebulizer chamber into very fine mist. The mains cable (6) is plugged into the device socket (14). Ensure that the values comply with those on the notice sign. The ultrasonic nebulizer is switched on with the main switch ON/OFF switch (13) at the rear of the device. The aspirated air is transferred to the nebulizer chamber (2) through the bacterial filter (5) via the 300 mm hose. The mist volume of the fluid generated in the nebulizer chamber (2) by means of the quartz oscillator is transferred near the patient's breathing organs by means of the generated air flow using the patient hose (4).

ON/OFF push button:


- Switch the device on using the on/off switch [1] on the front of the device.
- Green LED [2] illuminated/field label: I/O
- The device starts nebulizing

The red LED [3] lights up in case of faults

Mist volume push button:

Regulates the mist volume generated in the nebulizer chamber per unit of time.


Increasing the mist volume

- Press button (+) [1]
- Another green LED [2] lights up

Reducing the mist volume

- Press button (-) [1]
- A green LED [2] goes off

Air volume push button:


Increasing the air volume

- Press button (+) [1]
- Another green LED [2] lights up

Reducing the air volume

- Press button (-) [1]
- A green LED [2] goes off

Setting the timer:

Turning on the timer

- Set the desired timer time with the (+) or (-) buttons, see [2]
A maximum of 60 minutes can be set.
- If the timer is active, the green signal light is illuminated and displays the set time
When the set time has elapsed, the device switches off automatically:
- The red light [3, see note on ON/OFF push button] is lit and a signal tone can be heard. The signal tone stops after approx. 2 minutes.
- The device can only be switched back on after it has been switched off with the ON/OFF switch


Do not exceed the max. time of 15 minutes if sterile water capsules without a refilling system are used

Heatable patient hose

Only use original hoses. Malfunctions may arise if other heated hoses are used.

The hose heater switches on automatically when the ultrasonic nebulizer starts. The operating temperature of max. 40 °C is reached after approx. 10 minutes

- Connect the plug of the hose heater (E) to the socket at the rear of the device. Make sure the mark on the plug is pointing upwards.
- Attach the heated silicone hose (instead of the long 1000mm hose) onto the opening of the nebulizer chamber lid, the automatic feeder set or the medication container.

The connection pieces of the respective chamber lids have conical 22mm ISO standard connections, so that conventional hoses with these dimensions can be connected in addition to those listed in the appendix.

Caution, risk of patient injury

Only use heated hoses up to an ambient temperature of 30 °C

6 Putting into operation

6.1 Preparation

Place the device upright on a level, solid surface. The device must be upright to ensure proper function of the overflow protection. On a soft surface (e.g. pillow, mattress), the ventilation slots may be blocked, which can cause the device to overheat.

6.2 Connecting to power supply

Warning, risk of electric shock


Connecting devices to auxiliary power sockets may increase leakage current. If the protective ground on one of these devices fails, the leakage current may exceed the permissible values. Only make connections after agreement with the respective manufacturer of the device. Have service personnel check the leakage current. If the permissible value is exceeded, do not use the auxiliary power sockets on the device but separate mains power sockets on the wall.

Caution, device inoperable

Before putting the device into operation for the first time, the device voltage must be set to the mains voltage available using the voltage switch at the rear of the device. 230V is preset.


Selecting the mains voltage

- Carefully unlock the fuse insert (A) on the device socket using a screwdriver and remove it
- Remove the voltage switch (B) from the fuse insert and select the required mains voltage by rotating it. Put the fuse insert back in by pressing it until it audibly clicks into place.
- The required mains voltage (110 V/127 V/230 V) is now visible in the window (C)

Observe the rating plate on the underside of the device.

- Connect the power cable (D) to the device socket and to a grounded mains power socket.

6.3 Functional test

Check the functioning of the device prior to every use and after every reprocessing procedure.


Only device parts in faultless condition can ensure the functionality of the ultrasonic nebulizer. For this reason, the device parts must be subjected to a thorough check and inspection prior to any assembly.

Check and inspect:

- if the mains cable is undamaged,
- if the device parts made from plastic and rubber (e.g. hose, nebulizer chamber) are in faultless condition and do not show signs of age-related damage,
- if the device was properly cleaned and there is no residue or soiling.


6.4 Preparing USV 2010

6.4.1 Nebulizing directly from nebulizer chamber

Inserting the nebulizer chamber


- Fill the nebulizer chamber to the max. mark [2] with nebulizing fluid (2).
- Plug the nebulizing chamber (2) into the holder of the USV 2010 [1]. Observe the preferred position.
- Press the lid onto the nebulizer chamber

Mounting the bacterial filter

Connect the bacterial filter [1] to the fan connection [2] on the housing


Observe the flow direction according to the manufacturer specifications.

Mounting the 300 mm hose

A 300 mm hose (3) connects the bacterial filter (5) and the nebulizer chamber lid (18).

Connect the short hose [1] to the nozzle on the bacterial filter [2] and the nozzle on the lid [3].


Mounting the 1200 mm hose

Fit the long hose [1] to the second nozzle on the lid [2] and push it into the clips on the hinged arm. Guide the hose upwards in such a way that condensate can flow back into the nebulizer chamber.

Secure the lid of the nebulizer chamber with the holding-down device.

6.4.2 Nebulizing with standard disposable sterilized water systems

Fill the nebulizer chamber (2) with distilled water until the sterile water system [2] (e.g., Hudson RCI) is wetted up to the mark for contact water (shaded area [1]). The distilled water in the nebulizer chamber serves to transmit the oscillations to the sterile water system.

Prepare the sterile water system as specified in the respective instructions for use and fit it onto the nebulizer chamber in place of the lid.

Connect the hoses and secure the sterile water system in position with the holding-down device.

The contact water must be changed after every use, but at least once a day.

Sterile water capsules can be purchased from pharmacies or medical supplies retailers

Make sure that the nebulizer chamber and the sterile water system contain a sufficient amount of liquid. If this is not the case, this can damage the sterile water capsule (loss of sterility) or the nebulizer chamber. Always secure sterile water capsules with the holding-down device (29).

Warning, risk of injury!

If the sterile water capsule is empty, the water in the nebulizer chamber can heat up to more than 70 °C

6.4.2.1 Quartz module with thermal shutdown (recommended)

If sterile water capsules are used, the use of a quartz module with thermal shutdown is recommended for safety reasons. This can prevent the chamber and the sterile water capsule from being damaged. This also protects the safety of the patient and the user, as they cannot get injured by the heated water.

After the quartz module has cooled down, the device is ready for use after approx. 10 minutes.

6.4.3 Nebulization from the disposable medication container

- For nebulizing small quantities of medication.

Only use approved medication containers.


Use only disposable hoses and disposable medication containers for medication nebulization. Dispose of the hoses and the medication containers after single use.

Medication nebulization must only be carried out as directed by a medical professional.

Follow the instructions of the medication manufacturer regarding the suitability of the medication for nebulization.

Mounting the medication container

Special insert (Optional accessories) for the nebulizer chamber to nebulize medication

Fill the nebulizer chamber with distilled water until the sterile water system is wetted up to the mark [1] for contact water (shaded area). The distilled water in the nebulizer chamber serves to transmit the oscillations to the medication container.

Fill the medication into the medication container [2], cover it with the lid [3] and place it in the nebulizer chamber.

Connect the hoses and secure the medication container in position with the holding-down device.

The contact water must be changed after every use, but at least once a day.


Make sure that the nebulizer chamber contains a sufficient amount of liquid.

Do not allow the medication container to run dry

6.4.4 Nebulization with refilling system (optional)

Refilling systems permanently supply liquid to fill the nebulizer chamber.

Place the prepared refilling system [3] onto the nebulizer chamber instead of the lid. Connect the hoses and secure the refilling system in place with the holding-down device.

The height difference between the nebulizer chamber and the container (bag/bottle) must not exceed 1000 mm.

Pierce the container with the transfer set, e.g., Intrafix Air by B. Braun. Do not kink the hose.

The container must be ventilated maintaining sterile conditions.

Functional test of refilling system

When the container has been connected, check operation of the float:

- Fluid runs into the nebulizer chamber.
- Feed stops automatically when the fluid has reached its level (range between min. and max. mark).

In the event of a fault, clean the float valve and check operation again. Then put the ultrasonic nebulizer into operation.

End of operation of refilling system

After operation, lock the feed of the container, e.g., by closing the knurled screw on the transfer set.


6.5 Assembling accessories

Installation to the basic device must be in accordance with the instructions for use of the basic device on which this medical device is used. Make sure that the connection to the basic device system is secure.

Strictly observe assembly instructions and instructions for use.

6.5.1 Mounting the table model (Optional accessories)

- Fit the pole (300 mm long) into the holder at the back of the nebulizer and secure it in place with the locking screw.
- Slide the holding-down device (to secure the lid of the nebulizer chamber) and the hinged arm onto the pole and position them as required.

Avoid a sagging of the patient hose, see 6.4



Note

Now place the adjusting ring of the bottle holder (28) onto the stainless steel pole and fix it in place in the desired position using the rotary knob. Hang a container with the liquid to be nebulized on the bottle holder.

The air intake opening of the container should be opened so that the infusion can continue to flow into the nebulizer chamber.

6.5.2 Mounting the wall-mounted model accessory set (Optional accessories)



- 26 Stainless steel pole 25 x 600 mm
- 27 Hinged arm with adjusting ring
- 28 Bottle holder with adjusting ring
- 29 Holding-down device with adjusting ring
- 30 Rail clamp for standard rail

- Attach the bracket and the pole (600 mm long) to a standard rail (25 x 10 mm) on the wall.
- Fit the USV 2010 onto the pole using the holder at the back of the nebulizer and secure it in place with the locking screw.
- Slide the holding-down device (to secure the lid of the nebulizer chamber) and the hinged arm onto the pole and position them as required.


6.5.3 Mounting the trolley accessory set (Optional accessories)

- 31 Stainless steel pole
- 32 Hinged arm with adjusting ring
- 33 Bottle holder with adjusting ring
- 34 Holding-down device with adjusting ring
- 35 3-foot trolley

- Slide the retaining ring (A) onto the pole and screw tight.
- Fit the USV 2010 onto the pole of the trolley using the holder at the back of the nebulizer and secure it in place with the locking screw.
- Slide the holding-down device (to secure the lid of the nebulizer chamber), the hinged arm, and the bottle holder (for standard refilling systems) onto the pole and position them as required.

6.5.4 Transport with trolley


Warning, risk of tipping over!

For intrahospital transport push the ultrasonic nebulizer down to the pre-mounted stop (maximum 70 cm above the floor) and move it to the preferred position.
Guide the device during transport.

Caution

Risk of crushing during assembly, start-up, and use of the hinged arm.

When applying the lockable casters [B], make sure that the ultrasonic nebulizer is positioned between the lockable casters.

7 Operation

7.1 Functional test

Warning

The ultrasonic nebulizer is used for treating patients only in accordance with its intended use.
A performance restriction of the device can cause complications during the treatment.

The ultrasonic nebulizer must undergo a complete functional test before every use. Ensure that

- the bacterial filter is installed and functional,
- the hoses and the nebulizer chamber are installed and functional and that there is no residue,
- the hoses and the nebulizer chamber have a firm fit, seal tightly and that no mechanic forces are applied to the hoses; the hoses must not be kinked,
- device parts made from plastic and rubber are in faultless condition and do not show signs of aging,
- the mains cable is not damaged,
- the nebulizer chamber is filled sufficiently,
- the holding-down device rests on the respective lid.

7.2 Starting operation

Aerosol therapy must only be carried out as directed by a medical professional.

- Direct the unconnected end of the long hose to the mouth/nose area of the patient or to the opening of the tracheal cannula
- Maintain a minimum distance of 20 cm to the patient
- Set the toggle switch at the rear to "I" (on)
- Set the required nebulizer output
- Set the required air flow
- Make sure that there is enough fluid


Do not allow the nebulizer chamber to run dry

If the device runs dry, it switches off automatically and a signal tone can be heard. The signal tone stops after approx. 2 minutes.

If the red light indicates a malfunction:

- Turn the device off with the ON/OFF push button
- Set the toggle switch at the rear to "0" (off). The red light goes out
- For troubleshooting see 8 Faults, troubleshooting on page 29.
- Set the toggle switch back to "I" (on) and switch on the device with the ON/OFF push button.

7.3 Stopping operation

- Press push button [1]
- All LEDs are off and the device stops all functions
- Set the toggle switch at the rear to "0" (off)
- Close the hose clip of the refilling system
- If not in use (for more than 60 min.) Disconnect the mains plug

8 Faults, troubleshooting

8.1 General

The ultrasonic nebulizer underwent a thorough quality check at the factory. If a fault occurs despite these tests, you may be able to rectify them yourself if you observe the following notices.

There is a fault indicator on the left side of the control panel at the ON/OFF switch. This lights up red if a fault occurs or the set time has elapsed (timer end).

After completing the set time see 5.6

After a fault occurred, switch the device off.

Press the ON/OFF switch to do this. Remove the cause of the fault if possible and switch the device back on. The following can be potential causes of the fault:

Fault	Cause	Remedy
Red light is lit		
Not possible to switch on the device	The plug is not positioned correctly in the device or power socket	Check position of the plug in the device or power socket
	On/off switch at the rear of the device is not set to On (I)	Set on/off switch to On (I)
	On/off switch is switched off	Set switch to On
	No, defective or wrong mains fuse	Check main fuse Replace fuse as described in 8.3
	Electronics faulty	Have it repaired
	Quartz module with thermal shutdown	Thermal switch has triggered. The device is functional again after cooling down for 10 minutes
Device is not working/ Red light is lit	Nebulizer chamber not fitted correctly	Insert nebulizer chamber correctly. Observe the preferred position
	Nebulizer chamber empty	Refill nebulizer fluid
	Leak in nebulizer chamber	Screw quartz oscillator firmly into the nebulizer chamber
	Faulty quartz oscillator	Replacing the quartz oscillator
	Quartz module with thermal shutdown	Thermal switch has triggered. The device is functional again after cooling down for 10 minutes
	Electronics faulty	Have it repaired
Not enough nebulization when set at maximum output	Nebulizer fluid does not match scale	Refill nebulizer fluid or remove some
	Resistance of filter too high	Replace filter
	Condensate in hose	Drain condensate. Re-route hose in such a way that condensate can flow back into the nebulizer chamber.
	Leak in hoses or in the lid of the nebulizer chamber	Fit lid tightly and secure with holding-down device. Connect hoses firmly to nozzles.
Aerosol heating does not work	Plug of hose heater not connected	Connect plug
	Faulty hose	Replace hose
	Electronics faulty	Have it repaired
	Switched on for less than 10 min.	Reaching the operating temperature can take up to 10 min
No nebulization	Quartz oscillator failure	Replace quartz module
	No contact water	Refill contact water. Check chamber for leaks
	Accumulated condensation in the nebulizer hose prevents aerosol flow	Drain condensate. Perform installation as described in 6.4
	Bacterial filter or hoses not installed correctly	Perform installation as described in 6.4

8.2 Replacing the quartz oscillator

Screw the new quartz module (C) into the nebulizer chamber and check it for tight fit. The edge of the nebulizer chamber must be flush with the quartz module.

8.3 Replacing the fuses

Replace the fuses in the device plug as follows:

Warning, electric shock!

Disconnect the device from the mains before replacing the fuses.


1. Unplug the mains cable from the device socket (1).
2. Carefully unlock the lateral locks (3) on the device socket (1) using a small screwdriver and remove the fuse insert (4).
3. Remove defective fuses (2) and replace them with new ones.
4. Insert the fuse insert (4) and press it in until it audibly clicks into place.
5. Check tight fit of the fuse insert (4).

- | | |
|---|---------------|
| 1 | Device socket |
| 2 | Fuses |
| 3 | Locks |
| 4 | Fuse insert |

Warning, observe fuse type!

Only fuses of the following type may be used:
IEC 60127-2/ III T1L 250V

9 Reprocessing

Warning, electric shock!

Disconnect the device from the mains before replacing the fuses.

9.1 Disassembly

- Disconnect the mains plug and the device plug
- Disconnect the plug of the hose heater
- Disconnect the hoses and the bacterial filter.
- Push the holding-down device up and remove the lid of the nebulizer chamber. Detach the nebulizer chamber and empty out the fluid. Observe the hospital hygiene regulations
- Dispose of disposable articles (disposable hoses, medication container and lid, containers of sterile water systems) in an environmentally-friendly manner. Observe the hygiene regulations.
- Do not disinfect or sterilize the filter
- Replace and dispose of the filter as directed by the manufacturer.

9.2 Notes on implementation

The individual parts of the ultrasonic nebulizer can be cleaned and disinfected according to the information in the table below. Please observe the following information to do this:

- The parts must be cleaned and dried before every disinfection.
- Dirt particles can be encapsulated and the device parts will therefore not be sterile after disinfection. This is why all device parts must be thoroughly cleaned and persistent soiling removed with a cloth and cleaning agent.
- Please observe the instructions of the manufacturer and the hygiene expert before using disinfectants.

Warning, risk of injury!

Disinfectants can contain substances that are harmful to health and cause injuries when they come into contact with skin or eyes. Protect your skin and eyes and comply with the hygiene rules when working with disinfectants.

Caution, material changes!

Almost all components of the ultrasonic nebulizer are made of plastics. Solvents, some disinfectants and cleaning agents can corrode plastics or cause stress cracks. Do not use agents containing hydrocarbons or alcohol to clean the surfaces.

Caution, functional inspection!

The functionality of the device parts must be checked after every cleaning or disinfection (7.1 Functional test).

Note

If wrapping cloths that are not color fast are used, this can discolor the plastic parts.

Note

After every cleaning, disinfection and complete drying, check the rubber seal of the quartz module thoroughly for potential damages.

9.3 Implementation

Components	Cleaning	Disinfection ²	Dishwasher	Hot steam sterilisation ³
Nebulizer chamber ¹ (brown color)	in warm water with mild household cleaner	immerse in diluted solution	rinse cycle up to 95 °C	up to 134 °C 18min
Nebulizer chamber lid ¹ (white)				
Nebulizer chamber ¹ with quartz module				
Hoses			not allowed	
Nebulizer chamber ¹ (transparent)	in warm water with mild household cleaner	immerse in diluted solution	rinse cycle up to 95 °C	up to 121 °C 20 min
Nebulizer chamber lid ¹ (transparent)				
Bacteria and virus filter	Disposable items: replace after no later than 1 week or in accordance with particular manufacturer instructions. If a closed sterile water system is used, it is recommended to replace the filter after 48 hours of operation or every 14 days.			
Medication container	intended for single use			

Medication container lid			
Housing			
Power cable	wipe with moist cloth	wipe disinfection	not allowed
Accessory sets			
<ol style="list-style-type: none"> 1. Rinse the nebulizer chamber and nebulizer chamber lid with distilled water after every use. Disinfection must be performed with common medical disinfectants. The quartz module remains tightly screwed into the nebulizer chamber. 2. After the specified contact time and with the solution diluted (according to manufacturer's specifications), the components must be cleaned from residual disinfectant with a moist cloth and then dried. 3. Caution, hot steam accelerates the natural aging of plastics. The material change can affect the functionality of the device parts. 			

Examples of surface disinfectant:

- Incidin® Plus *
- Incidin® Perfect *
- Antiseptica Combi Surface *
- Antifect® FF *
- B 10 *

* registered trademarks of the respective manufacturers

Examples of device disinfectants:

- Lysetol® AF *
- Gigasept® FF *
- Sekusept® forte S*

All individual parts must be rinsed with distilled water after every use.

Disposable items must not be reprocessed due to risk of infection. The above-mentioned examples for disinfectants are suggested by the device manufacturer. The recommendation of the listed active ingredient basis (pH 4.5-8) applies. This does not release the user from consulting the hygiene expert in individual cases or checking the components to be disinfected for compatibility with the active ingredients.

9.4 Manual cleaning of the quartz

For cleaning, remove the nebulizer chamber lid (18) from the nebulizer chamber (17) and screw the quartz module (20) out of the nebulizer chamber (17). Reassembly is performed in reverse order.

9.5 Cleaning the automatic feeder set

Information can be found in the separate instructions for use

10 Preventive maintenance / repair**10.1 General****Caution, risk to health!**

The ultrasonic nebulizer is used for treating patients. The device or parts thereof can be contaminated. This is why, prior to returning it for inspection or repair, the nebulizer chamber, the automatic feeder set, the bacterial filter and all hoses must be removed and the device must be cleaned and disinfected.

10.2 Preventive maintenance

The ultrasonic nebulizer is maintenance-free. The operating personnel is not required to perform preventive maintenance other than normal care. To ensure functionality of the product, visual and functional tests must be performed prior to use. In order to ensure operational safety of the product and availability of all functions as well as increased life span, we recommend yearly preventive maintenance performed by an authorized specialist dealer/service technician. The scope of any preventive maintenance depends on a device-specific performance plan. For this, a safety check in accordance with the standard IEC 601-1 is performed.

The distributor prescribes a yearly safety check in accordance with §6 MPBetreibV (German Medical Device Operator Ordinance).

10.3 Repair**Disclaimer**

If unauthorized persons interfere with the device, no warranty claims can be made. The repair may only be performed by service technicians authorized by the distributor. These authorized service technicians can request wiring diagrams, descriptions, spare parts lists and test plans.

If they detect any faults, the product must no longer be used. In this case, record the serial number on the rating plate and contact the service address (see cover page) or your authorized specialist dealer.

11 Available products**11.1 Standard delivery**

consisting of:

2010-90-0-000

- Ultrasonic nebulizer, basic device
- Power cable
- Sterilizable nebulizer chamber, complete with nebulizer chamber lid and quartz module
- Disposable corrugated hose approx. 300 mm
- Disposable corrugated hose approx. 1150mm
- Bacterial filter

11.2 Optional accessories

The following accessories are not included in the scope of delivery and must be ordered separately:

Trolley accessory set

2002-00-5-900

consisting of:

- 5-foot trolley, chrome-plated, 50 mm wheels,
- Stainless steel pole, hinged arm with adjusting ring,
- Bottle holder with adjusting ring (stainless steel),
- Holding-down device with adjusting ring (stainless steel)

Table model accessory set

2002-01-3-000

consisting of:

- Stainless steel pole 25 x 600 mm, hinged arm with
- adjusting ring, bottle holder with adjusting ring (stainless steel),
- Holding-down device with adjusting ring (stainless steel)

Wall-mounted model accessory set

2002-01-4-000-00

consisting of:

- Stainless steel pole 25 x 600 mm, hinged arm with
- adjusting ring, bottle holder with adjusting ring (stainless steel),
- Holding-down device with adjusting ring (stainless steel), rail clamp for standard rail 25 x 10 mm (DIN)

Patient hose (silicone hose), heatable

2100-10-0-000

Length approx. 1100 mm, with circular connector, without holes for temperature sensors

11.3 Spare parts and consumables

Bacterial filter, complete, disposable	2150-00-0-001-00
Medication container, container and lid	2002-02-2-900
Corrugated hose 300 mm	2002-02-1-928
Corrugated hose 1200 mm	2002-02-1-932
Nebulizer chamber, complete	2002-02-1-684
Nebulizer chamber lid	2002-00-1-902
Nebulizer chamber	2002-02-1-952
Quartz module with thermal switch	2002-02-1-912
Hinged arm	5003-24-0-000
Holding-down device, stainless steel	3151-02-1-000
Bottle holder, stainless steel	3101-02-1-000

**11.4 Interface description**

All devices or accessories combined with the ultrasonic nebulizer must be listed in 11.2 and 11.3 and meet the requirements of the following interface description in accordance with Zweckbestimmung under the intended operating conditions. Only accessories and consumables suitable for the specific application, which have with CE symbol and possess the characteristics defined in the interface description, are approved. The authorized personnel is responsible for configuration of the overall system and functionality checks.

12 Disposal

At the end of its life span, the ultrasonic nebulizer must be properly disposed of. Take care to separate the material accurately. The device and accessories do not contain any hazardous materials; the device component materials can be fully recycled. This product falls under the scope of EC Directive 2002/96/EC (Waste Electrical and Electronic Equipment Directive). The distributor or an authorized dealer takes back the device at the end of its use period for a fee. Disposal via the municipal collection points for waste electrical and electronic equipment is not allowed. For more information on legal disposal, please contact your authorized specialist dealer.

Outside of the European Economic Community

For disposal of this product, the applicable national regulations regarding disposal and handling of waste electronic equipment must be complied with.

13 Technical data**13.1 Ambient conditions**During operation:

Temperature	10 °C to 35°C
Temperature when using the hose heater with heatable hose	10 °C to 30 °C
Relative humidity	0 to 90 % relative humidity, non-condensing

During storage:

Temperature	-20 °C to 70 °C
Relative humidity	0 to 90 % relative humidity, non-condensing

13.2 Performance data

Voltage	110V/ 127V/ 230V
	50 - 60Hz
Fuses	IEC 60127-2/ III T1L 250V, 2 pieces
Power consumption:	
Ultrasonic nebulizer	50VA
Hose heater	30VA
Resistance of hose heater	3 Ohm ± 10 %
Ultrasonic frequency	1.68 MHz ± 5 %

Sound pressure level	max. 35 dB(A)
Protection class	Protection class I, type B
Protection class	IP 21 (dripping water)
Operating mode	Continuous operation, time-controlled
Nebulizer output	max. 3mL/min
Nebulizer output with sterile water system	up to 2.5 mL/min
Nebulizer chamber volume	0 – "Min" approx. 80 mL; "Min" – "Max" approx. 85 mL; Total volume approx. 170 mL
Size of droplets	86 % < 5 µm, average 4 µm
Air flow	up to 20 L/min if the bacterial filter is used
Device dimensions	
W x H x D	205 x 190 x 315 mm
Weight	approx. 3.5 kg
Classification according to EC-Directive 93/42/EEC	Class IIa
Annex IX	
UMDNS Code	12-719

14 Electromagnetic compatibility

14.1 Electromagnetic emission

Guideline and manufacturer's declaration – electromagnetic emissions		
The ultrasonic nebulizer is intended for use in an electromagnetic environment as specified below. The customer or user of the product should ensure that it is operated in such an environment.		
Interference	Corresponds to	Electromagnetic environment – guideline
RF emissions according to CISPR11	Group 2	The ultrasonic nebulizer generates RF energy internally and transfers this to the transducer for mist generation. The RF emission is very low. However, it can cause interference in nearby electronic equipment.
RF emissions according to CISPR11	Class B	The ultrasonic nebulizer is suitable for use in all clinical facilities, including those in residential areas (specialist personnel), and those directly connected to a public supply grid that also supplies buildings.
Harmonic emissions in accordance with IEC 61000-3-2	Class B	
Voltage fluctuations/flicker emissions in accordance with IEC 61000-3-3	corresponds	

14.2 Electromagnetic interferences in other devices

Guideline for avoiding, detecting and remedying electromagnetic interferences in other devices		
Other electrical/electronic devices should not be operated near the ultrasonic nebulizer or stacked on top or below it. If such an arrangement cannot be avoided, the respective devices should be monitored to check proper function. Due to the various device functions, the effects can differ significantly and may be hard to detect. Examples:		
Device	Fault	Remedy
Radio and TV devices	- noise / crackling sound - horizontal stripes on the screen	- increase distance
Monitoring systems e.g. baby monitors	- noise / crackling sound	- change arrangement
Cordless phones	- noise / crackling sound	- change orientation
Wireless thermometer Wireless weather stations	- interference during data transmission - missing or wrong display	- switch receiver channel
General electronic devices	Malfunction e.g. stop or change of the intended operation mode	

14.3 Electromagnetic immunity

Guidelines and manufacturer's declarations – electromagnetic immunity			
The ultrasonic nebulizer is intended for use in the electromagnetic environment specified below. The customer or user of the product should ensure that it is operated in such an environment.			
Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment – guideline
Electrostatic discharge/ ESD according to IEC 61000-4-2	± 6 KV contact discharge ± 8 KV air discharge	± 6 KV contact discharge ± 8 KV air discharge	Floors should be made from wood or concrete or tiled with ceramic tiles. If the floor is covered with synthetic material, the relative humidity level must be at least 30 %.
Electrical fast transients/bursts according to IEC 61000-4-5	± 2 KV for power supply lines	± 2 KV for power supply lines	Mains voltage quality should be that of a typical commercial or hospital environment.

Surges on AC mains lines according to IEC 61000-4-11	± 1 KV differential mode voltages ± 2 KV common mode voltages	± 1 KV differential mode voltages ± 2 KV common mode voltages	Mains voltage quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and fluctuations of the supply voltage according to IEC 61000-4-11	>95 % dip, 0.5 periods 60 % dip, 5 periods 30 % dip, 25 periods >95 % dip, 5 seconds	-100 % -60 % -30 % -100 %	Mains voltage quality should be that of a typical commercial or hospital environment. If the user of the ultrasonic nebulizer requires continued operation even during energy supply interruptions, it is recommended that the ultrasonic nebulizer is powered from an uninterruptible power supply or a battery.
Power frequency magnetic field (50/60 Hz) according to IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical commercial or hospital environment.
Note: U_T is the alternating mains voltage prior to application of test levels.			

14.4 Electromagnetic immunity, devices not intended for life support

Guidelines and manufacturer's declarations – electromagnetic immunity, devices not intended for life support			
The ultrasonic nebulizer is intended for use in the electromagnetic environment specified below. The customer or user of the product should ensure that it is operated in such an environment.			
Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment – guideline
Conducted RF disturbances according to IEC 61000-4-6	3 V _{eff} 150 kHz – 80 MHz	3 V _{eff} 150 kHz – 80 MHz	Portable and mobile radio equipment should not be used at a smaller distance to the product, including lines, than the recommended separation distance, which is calculated by means of the equation appropriate for the transmission frequency. Recommended separation distance: $d = 1.2 \sqrt{P}$ (for 150 kHz – 80 MHz) $d = 1.2 \sqrt{P}$ (for 80MHz – 800 MHz) $d = 2.4 \sqrt{P}$ (for 800 MHz – 2.5 GHz) with P as nominal output of the transmitter in watts (W) in accordance with the transmitter manufacturer's specifications and d as recommended separation distance in meters (m) In accordance with an on-site examination (see a), the field strength of fixed radio frequency transmitters should be lower than the compliance level (see b) for all frequencies. Interferences are possible near devices with the following symbol. 
Radiated RF disturbances according to IEC 61000-4-3	3 V/m 80 MHz – 2.5 GHz	3 V/m 80 MHz – 2.5 GHz	
Note 1: At 80 MHz and 800 MHz the higher frequency range applies Note 2: These guidelines may not be applicable in all cases. The propagation of electromagnetic variables is influenced by absorptions and reflections of the building, objects and persons. a) The field strength of fixed transmitters, such as base stations of cordless phones and mobile terrestrial radio equipment, amateur radio stations, AM and FM radio and TV channels, can in theory not be accurately predetermined. In order to assess the electromagnetic environment regarding the fixed transmitters, an electromagnetic site survey should be considered. If the measured field strength at the location where the ultrasonic nebulizer is operated exceeds the above-mentioned compliance levels, the device should be monitored to verify normal operation. If abnormal performance features are observed, additional measures may be required, e.g., a different location of the ultrasonic nebulizer. b) For the frequency range of 150 kHz to 80 MHz, the field strength should be lower than 3 V/m.			

14.5 Recommended separation distances

Recommended separation distances between portable and mobile RF telecommunication equipment and the product
The ultrasonic nebulizer is intended for operation in an electromagnetic environment, in which the RF disturbances are controlled. The customer or user of the product can help avoid electromagnetic interferences by complying with the minimum distance between portable and mobile RF telecommunication equipment (transmitters) and the product –depending on the output power of the communication device, as detailed below.
Separation distance in meters depending on the transmission frequency

Nominal output of the transmitter in watts	150 KHz – 80 MHz	80 MHz – 800 MHz	800 MHz – 2.5 GHz
	$d = 1.2 \sqrt{P}$	$d = 1.2 \sqrt{P}$	$d = 2.4 \sqrt{P}$
0.01	0.12	0.12	0.24
0.10	0.38	0.38	0.76
1.00	1.20	1.20	2.40
10.00	3.80	3.80	7.60
100.00	12.00	12.00	24.00

For transmitters whose maximum nominal output is not listed in the above table, the recommended separation distance d in meters (m) can be determined using the equation of the respective column, where P is the maximum nominal output of the transmitter in watts (W) according to the transmitter manufacturer's specifications.

Note 1:

At 80 MHz and 800 MHz the higher frequency range applies

Note 2:

These guidelines may not be applicable in all cases. The propagation of electromagnetic variables is influenced by absorptions and reflections of the building, objects and persons.

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© **All rights reserved**
by W. Krömker GmbH
subject to change without notice
Jun. 2016

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von W. Krömker GmbH
Drägerwerk AG & Co.KGaA